SUPABON SUTURE

STERILE SYNTHETIC ABSORBABLE SUTURE

Description :

SUPABON suture is a synthetic braided absorbable sterile surgical suture composed of a braided polyglycolic acid polymer . SUPABON is available in a range of gauge sizes and lengths attached to stainless steel needle of varying types and sizes. Full details are contained in the catalogue . Color is violet but seldom is green.

Indications :

SUPABON sutures are intended for use in general soft tissue approximation and / or ligation, including use in ophthalmic surgery , peripheral nerve anastomosis and microsurgery for vessels less than 2 mm in diameter . The safety and effectiveness of SUPABON sutures in cardiovascular tissue have not been established .

Application:

Sutures should be selected and implanted depending on patient condition , surgical experience , surgical technique and wound size.

Performance:

SUPABON suture elicits a minimal initial inflammatory reaction in tissues and ingrowth of fibrous connective tissue . Progressive loss of tensile strength and eventual absorption of SUPABON sutures occurs by means of hydrolysis. Absorption begins as a loss of tensile strength followed by a loss of mass. All of the original tensile strength is lost between four and five weeks post implantation . Absorption of SUPABON suture is essentially complete between 56and 70 days.

Days after Implantation	Approximate % original Strength Remaining
14 days	75%
21days(6-0and larger)	50%
21 days (7-0and smaller)	40%

Contraindications:

These sutures , being absorbable should not be used where extended approximation of tissues under stress is required .

Warnings / Precautions / Interactions:

Users should be familiar with surgical procedures and techniques involving absorbable sutures before employing SUPABON suture for wound closure , as risk of wound dehiscence may vary with the site of application and the suture material used. Surgeons should consider the in vivo performance (under PERFORMANCE section)when selecting a suture. As with any foreign body , prolonged contact of any suture with salt solutions, such as those found in the urinary or biliary tracts , may result in calculus formation . As an absorbable suture SUPABON may act transiently as a foreign body . Acceptable surgical practice should be followed for the management of contaminated or infected wounds . As this is an absorbable suture material , the use of supplemental non absorbable sutures should be considered by the surgeon in the closure of the sites which may undergo expansion , stretching or distension , or which may require additional support .

Skin sutures which must remain in place longer than 7 days may cause localized irritation and should be snipped off or removed as indicated. Under some circumstances, notably orthopedic procedures , immobilization of joints by external support may be employed at the discretion of the surgeon. Consideration should be taken in the use of absorbable sutures in tissues with poor blood supply as suture extrusion and delayed absorption may occur. Subcuticular sutures should be placed as deeply as possible to minimize the erythema and induration normally associated with the absorption process. This suture may be inappropriate in elderly , malnourished or debilitated patients , or in patients suffering from conditions which may delay wound healing . When handling this or any other suture material , care should be taken to avoid damage. Avoid crushing or crimping damage due to application of surgical instruments such as forceps or needle holders. Care should be taken to avoid damage when handling surgical needles . Garsp the needle in an area one –third (1/3) to one – half (1/2) of the distance from the attachment end to the point . Grasping in the point area could impair the penetration performance and cause fracture of the needle. Grasping at the butt or attachment end could cause bending or breakage. Reshaping needles may cause them to lose strength and be less resistant to bending and breaking . Users should exercise caution when handling surgical needles to avoid inadvertent needle stick injury, Discard used needles in Sharps containers . Adequate knot security requires the standard surgical technique of flat and square ties with additional throws as indicated by surgical circumstances and the experience of the surgeon . The use of additional throws may be particularly appropriate when knotting any monofilament suture.

Adverse Reactions:

Adverse reactions associated with the use of this device include transitory local irritation at the wound site , transitory inflammatory foreign body response, erythema and induration during the absorption process of subcuticular sutures . Like all foreign bodies SUPABON may potentiate an existing infection.

Sterility:

SUPABON sutures are sterilized by ethylene oxide gas . Do not resterilize!

Storage:

Recommended storage conditions:keep away from of the sunlight & temprature. Do not use after expiry date !

*Do not use if package is damaged.

*Destroy not used suture .

* Color chemical name : VIOLET: 1-hydroxy-4-[(4-methyl-phenyl)amino]-9,10-anthracenedione
GREEN: 1,4-bis[(methyl-phenyl)amino]-9,10-anthracenedione

Warning : If you re-use it will be in microbe pollution danger.

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