SURGICAL CHROMIC - CATGUT

STERILE ABSORBABLE

SUTURE (hardened strands)

Description :

Surgical catgut suture is a sterile absorbable surgical suture composed of purified connective tissue (collagen) derived from either the serosal layer of cattle (bovine) or the sub mucosal layer of sheep (ovine) intestines . Chromic catgut is surgical catgut which has been treated to prolong its resistance to digestion by adding the chromic salt solution . It is brown colored . Chromic catgut is available in a range of gauge sizes and lengths , attached to stainless steel needles of varying types and sizes . Chromic catgut is packed in conditioning fluid .Full details are contained in the catalogue.

Indications :

Chromic catgut is for use in general soft tissue approximation and/ or ligation , including use in ophthalmic procedures .

Application:

Sutures should be selected and implanted depending on patient condition , surgical experience , surgical technique and wound size.

Performance:

When chromic catgut suture is implanted , a moderate tissue inflammation occurs which is characteristic of a foreign body response . Loss of tensile strength and a loss of suture mass follows as the proteolytic enzymatic digestive process absorbs the surgical catgut . Due to the inherent variability of a natural material , figures for strength loss and absorption following implantation are for guidance only . Effective strength is usually lost in the range of 14 to 28 days and absorption is essentially complete by 90 days. In humans , various factors may affect the tensile strength loss and rate of absorption . Two of the major factors are :

Infection :

Chromic catgut is absorbed more rapidly in infected tissue than in non – infected tissue.

Tissue sites :

Chromic catgut is absorbed more rapidly in tissues where increased levels of proteolytic enzymes are present, as in the secretions from the stomach, cervix and vagina.

Contra- Indication:

These sutures , being absorbable , should not be used where extended approximation of tissues under stress is required . Chromic catgut is contraindicated in patients with known sensitivities or allergies to collagen since catgut is a collagen based material.

Warnings/Precautions / Interaction :

Users should be familiar with handling and knotting techniques involving catgut sutures before using this material for wound closure since the risk of wound dehiscence may vary with the site of application and the suture material used . Surgeons should consider the in-vivo performance of catgut (under performance section) when selecting this suture. This suture may be inappropriate in elderly , malnourished or debilitated patients or in patients suffering from conditions which may delay wound healing . As with any foreign body , prolonged contact of any suture with salt solutions , such as those found in the urinary or biliary tracts , may result in calculus formation . As an absorbable suture, chromic catgut acts transiently as a foreign body . Acceptable surgical practice should be followed for the management of contaminated or infected wounds. As with all absorbable suture materials, the use of supplemental non – absorbable sutures should be considered by the surgeon in the closure of sites which may undergo expansion , stretching , or distension or which may require additional support . Some patients may be hypersensitive to collagen and may exhibit an immunological reaction . Under some circumstances , notably orthopaedic procedures , immobilization of joints by external support may be employed at the discretion of the surgeon . Consideration should be taken in the use of absorbable sutures in tissues with poor blood supply as suture extrusion and delayed absorption may occur. In handling catgut or any other suture material , care should be taken to avoid damage . Avoid crushing or crimping damage due to application of surgical instruments such as forceps or needle holders. Adequate knot security requires the standard surgical technique of flat and square ties with additional throws as indicated by surgical circumstances and the experience of the surgeon. Avoid unnecessary tension when running down knots, to reduce the occurrence of surface fraying and weakening of the strand. Care should be taken when opening the pack as it may contain conditioning fluid. Care should be taken to avoid damage when handling surgical needles. Grasp the needle in an area one third (1/3) to one half (1/2) of the distance from the attachment end to the point. Grasping in the point area could impair the penetration performance and cause fracture of the needle . Grasping at the butt or attachment end could cause bending or breakage. Re- shaping needles may cause them to lose strength and be less resistant to bending and breaking . Users should exercise caution when handling surgical needles to avoid inadvertent needle stick injury . Discard used needles in Sharps containers.

Adverse Reactions :

Adverse reactions associated with the use of this device include : allergic response in certain patients and transient local irritation at the wound site followed by moderate transient inflammatory foreign body response . Like all foreign bodies catgut may potentiate an existing infection.

Sterility :

Chromic catgut sutures are sterilized by irradiation . Do not resterilize .

Storage:

Recommended storage conditions. keep away from direct sunlight & heat . Do not use after expiry date .

*Do not use if package is damaged.

*Destroy not used suture .

*Add. Material :

Tri ethanolamine	0.15%
Ethanol	90 %
Glycerol	2%

Warning : If you re-use it will be in microbe pollution danger.

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